

## qPCR 96-WELL PLATE COMPATIBILITY CHART

	MB-P96	MB-Q96	MB-Q96-S	MB-Q96-LP	MB-Q96-LBR	MB-Q96-LBR-W	Note / Holder NO.
Standard Well (0.2 ml, height 21mm.)	●	●	●				
Low Profile Well (0.1 ml, height 16mm.)				●	●	●	
<b>Applied Biosystems (Life Technology)</b>							
7000		●					
7300		●					4367032
7500		●					4367032
7500Fast				●			4359652
7900HT Fast Model(0.1 ml)				●			
7900HT Standard (0.2 ml)		●					
StepOne Plus				●			
ViiA7 Fast Model (0.1 ml)				●			4459846
ViiA7 (0.2 ml)		●					4459845
QuantStudio Fast Model (0.1 ml)				●			4459846
QuantStudio (0.2 ml)		●					4459845
QuantStudio 3/5 Fast (0.1 ml)				●			
QuantStudio 3/5 (0.2 ml)		●					
5700/7700		●	●				Sequencer
3100/3130 (xL) (0.1 ml Fast 96 Well)				●			4367470 (Retainer) 4315933 (Septa)
310/3100/3130 (xL) (0.2 ml Standard 96 Well)		●	●				4317241 (Retainer) 4315933 (Septa)
3700/3730 (xL) (0.1 ml Fast 96 Well)				●			4367469 (Retainer) 4315933 (Septa)
3700/3730 (xL) (0.2 ml Standard 96 Well)		●	●				4410228 (Retainer) 4315933 (Septa)
3500(xL) (0.1 ml Fast 96 Well)				●			4409530 (Retainer) 4412614 (Septa)
3500(xL) (0.2 ml Standard 96 Well)		●	●				4410228 (Retainer) 4412614 (Septa)
<b>Roche</b>							
LightCycler 480 & 480II						●	
LightCycler 96						●	
<b>Bio-Rad</b>							
iCycler iQ	●		●				
iQ5	●		●				
MyiQ	●		●				
Chromo4	●		●				
Opticon I & II					●	●	
CFX96					●	●	
<b>Stratagene/Eppendorf/Biometra</b>							
Mx4000/3000P/3005P	●		●				
AriaMx					●	●	
MasterCycle ep Realplex	●		●				
Toptical Thermocycler	●		●				